

Product Information	
Catalog No.	FC14923
Product name	Anti-Mouse CD360/IL21R Antibody (SAA0483), PE
Applications	ELISA, FCM
Species reactivity	Human, Mouse, Rat, Monkey
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	SAA0483
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL21R, NILR, Interleukin-21 receptor, IL-21R, IL-21 receptor, CD360, Novel interleukin receptor
Accession	Q9JHX3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.