

Product Information	
Catalog No.	FC15076
Product name	Anti-Human CADM3/NECL-1/TSLL1 Antibody (SAA1480), PE
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH-mFc
Clone ID	SAA1480
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	TSLL1, Immunoglobulin superfamily member 4B, Brain immunoglobulin receptor, NECL1, CADM3, Synaptic cell adhesion molecule 3, Cell adhesion molecule 3, IgSF4B, NECL-1, Nectin-like protein 1, TSLC1-like protein 1, SynCAM3, SYNCAM3, IGSF4B
Accession	Q8N126
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.