

Product Information	
Catalog No.	FC15078
Product name	Anti-Human HLA-B*27 Antibody (H4H2), PE
Applications	FCM
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	H4H2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	HLA-B, HLA class I histocompatibility antigen, B alpha chain, Human leukocyte antigen B, HLAB, HLA-B*27
Accession	P01889
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.