

Product Information	
Catalog No.	FC15093
Product name	Anti-Human Complement C5b-9 Complex Antibody (aE11#), PE
Applications	ELISA, FCM, WB
Species reactivity	Human, Baboon, Horse, Pig
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	aE11#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	C5b-9 Complex, Complement C5b-9 Complex, Complement Component 5b-9
Accession	C5b-9 Complex
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.