

Product Information	
Catalog No.	FC15100
Product name	Anti-Staphylococcus aureus Protein A/SPA Antibody (SAA1452), PE
Applications	ELISA, FCM, WB
Species reactivity	Staphylococcus aureus (P38507), Staphylococcus aureus (strain N315) (P99134), Staphylococcus aureus (strain Mu50 / ATCC 700699) (P0A015)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA1452
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	spa, Immunoglobulin G-binding protein A, Staphylococcal protein A
Accession	P38507
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.