

Product Information	
Catalog No.	FC15127
Product name	Anti-Human CD20/MS4A1 Antibody (SAA2209), FITC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA2209
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Purification	Affinity chromatography.
Target	B-lymphocyte surface antigen B1, Membrane-spanning 4-domains subfamily A member 1, MS4A1, Leukocyte surface antigen Leu-16, B-lymphocyte antigen CD20, Bp35, CD20
Accession	P11836
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.