

Product Information	
Catalog No.	FC15175
Product name	Anti-CD77 Antibody (38.13), APC
Applications	FCM
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	38.13
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Purification	Affinity chromatography.
Target	CD77, globotriaosylceramide (Gb3), Burkitt's lymphoma associated antigen (BLA), ceramide trihexoside (CTH), and Pk blood group antigen
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.