


Product Information	
Catalog No.	FC15218
Product name	Anti-TF-Ag Antibody (Karo11)
Applications	ELISA, FCM
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	Karo11
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	T-Ag, TF-Ag, T antigen, TF epitope, Thomsen-Friedenreich epitope, β -D-Galp-(1 $\rightarrow$ 3)-D-GalpNAc, Gal β 1-3GalNAc, (Gal)1 (GalNAc)1
Accession	CAS: 20972-29-6, CAS: 66798-32-5
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-TF-Ag Antibody (Karo11)

SDS-PAGE