

Product Information	
Catalog No.	MA10006
Product name	Anti-Human IgM Monoclonal Antibody (1A184), FITC
Applications	ELISA, FCM, IF, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A184
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Purification	Affinity chromatography.
Immunogen	Native Human IgM.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Alternative Names	IgM
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.