

Product Information	
Catalog No.	MA10016
Product name	Anti-ASFV K196R Antibody (1A009)
Applications	ELISA, WB
Species reactivity	African swine fever virus
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A009
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	Thymidine kinase, 2.7.1.21, K196R, K196R CDS, ASFV-Georgia_4-074
Immunogen	Recombinant ASFV K196R protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.