

Product Information	
Catalog No.	MA10033
Product name	Anti-ASFV B117L Antibody (1A026)
Applications	ELISA, WB
Species reactivity	African swine fever virus
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A026
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	hypothetical protein IM014_gp118, B117L CDS
Immunogen	Recombinant ASFV B117L protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.