


Product Information

Catalog No.	MA10059
Product name	Anti-Human B2M/ β 2-MG Antibody (1A053)
Applications	ELISA (Dec), GICA (Dec), CLIA (Dec)
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A053
Form	Liquid
Storage buffer	PBS buffer, pH 7.4, containing 0.1% NaN ₃ as preservative.
Purification	Protein A/G purified from cell culture supernatant.
Target	Beta-2-microglobulin, B2M
Immunogen	Human Beta-2-microglobulin Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human B2M/ β 2-MG Antibody (1A053)

SDS-PAGE