

Product Information	
Catalog No.	MA10062
Product name	Anti-Human CD284/TLR4 Monoclonal Antibody (1A057)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG3
Clone ID	1A057
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Toll-like receptor 4, TLR4, CD284, hToll,
Immunogen	Mammalian cells - derived recombinant Human CD284 / TLR4 (Met1-Lys631).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.