

Product Information	
Catalog No.	MA10081
Product name	Anti-Agrobacterium sp. aroA/EPSPS Monoclonal Antibody (1A076)
Applications	ELISA, IHC, WB
Species reactivity	Agrobacterium sp. (strain CP4)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A076
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	3-phosphoshikimate 1-carboxyvinyltransferase, 5-enolpyruvylshikimate-3-phosphate synthase, EPSP synthase, EPSPS, aroA, CP4-EPSPS, EPSPS-CP4
Immunogen	E. coli - derived recombinant Agrobacterium sp aroA / EPSPS (His3-Ala455).
Accession	Q9R4E4
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.