

Product Information	
Catalog No.	MA10085
Product name	Anti-Human GPC3 Monoclonal Antibody (1A080)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A080
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	GPC3, MXR7, GTR2-2, Glypican-3, Intestinal protein OCI-5, OCI5,
Immunogen	Mammalian cells - derived recombinant Human GPC3 (Met1-His559).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.