

Product Information	
Catalog No.	MA10090
Product name	Anti-Human WFDC2/WAP5 Monoclonal Antibody (1A085)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	1A085
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Epididymal secretory protein E4, Major epididymis-specific protein E4, HE4, WFDC2, WAP5, WAP four-disulfide core domain protein 2, Putative protease inhibitor WAP5,
Immunogen	Mammalian cells - derived recombinant Human WFDC2 / WAP5 (Met1-Phe124).
Accession	Q14508
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.