

| Product Information   |                                                                                                                                                               |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | MA10100                                                                                                                                                       |
| Product name          | Anti-Mouse IL6 Monoclonal Antibody (1A095)                                                                                                                    |
| Applications          | ELISA, IHC, WB                                                                                                                                                |
| Species reactivity    | Mouse                                                                                                                                                         |
| Host species          | Mouse                                                                                                                                                         |
| Clonality             | Monoclonal                                                                                                                                                    |
| Isotype               | IgG1                                                                                                                                                          |
| Clone ID              | 1A095                                                                                                                                                         |
| Form                  | Liquid                                                                                                                                                        |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                            |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                           |
| Target                | Hybridoma growth factor, BSF-2, IFN-beta-2, CDF, IL6, Interleukin-6, IL-6, B-cell stimulatory factor 2, Interferon beta-2, IFNB2, CTL differentiation factor, |
| Immunogen             | E. coli - derived recombinant Mouse IL6 (Phe25-Thr211).                                                                                                       |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.  |
| Note                  | For research use only.                                                                                                                                        |