

Product Information	
Catalog No.	MA10120
Product name	Anti-Human FLT1 Monoclonal Antibody (1A120)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A120
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	FLT1, Vascular endothelial growth factor receptor 1, VEGFR1, Vascular permeability factor receptor, FRT, Tyrosine-protein kinase receptor FLT, VEGFR-1, FLT, FLT-1, Tyrosine-protein kinase FRT, Fms-like tyrosine kinase 1
Immunogen	Mammalian cells - derived recombinant Human FLT1 (Met1-Ile328).
Accession	P17948
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.