

Product Information	
Catalog No.	MA10135
Product name	Anti-Human ATAT1 Monoclonal Antibody (1A143)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A143
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	C6orf134, TAT, MEC17, Alpha-TAT, ATAT1, Alpha-TAT1, Acetyltransferase mec-17 homolog, Alpha-tubulin N-acetyltransferase 1,
Immunogen	E. coli - derived recombinant Human ATAT1 (Met1-Pro196).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.