

Product Information	
Catalog No.	MA10140
Product name	Anti-Eimeria tenella Membrane-associated calcium-binding protein Monoclonal Antibody (1A148)
Applications	ELISA, IHC, WB
Species reactivity	Eimeria tenella
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A148
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Membrane-associated calcium-binding protein, related, putative, ETH_00009450
Immunogen	E. coli - derived recombinant Eimeria tenella Membrane-associated calcium-binding protein, related, putative (Gly24-Leu304).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.