

Product Information	
Catalog No.	MA10142
Product name	Anti-HSV2/HHV-2 gD/US6 Monoclonal Antibody (1A150)
Applications	ELISA, IHC, WB
Species reactivity	Human herpesvirus 2 (HSV-2/HHV-2)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A150
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	Envelope glycoprotein D, gD, US6
Immunogen	E. coli - derived recombinant HHV-2 gD / US6 (Lys26-Pro339).
Accession	P03172, Q69467
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.