

Product Information	
Catalog No.	MA10151
Product name	Anti-SARS-CoV-2 S2 Protein Monoclonal Antibody (1A160)
Applications	ELISA, IHC, WB
Species reactivity	SARS-CoV-2 (2019-nCoV)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	1A160
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	S2 Protein, Spike protein S2, Spike protein S2', Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1, S
Immunogen	E. coli - derived recombinant SARS-CoV-2 S2 Protein (Tyr1047-Gln1201).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.