

Product Information	
Catalog No.	MA10157
Product name	Anti-Human CD31/PECAM1 Monoclonal Antibody (1A168)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A168
Form	Liquid
Storage buffer	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium azide.
Purification	Protein A/G purified from cell culture supernatant.
Target	PECAM-1, Platelet endothelial cell adhesion molecule, GPIIA', PECAM1, PECA1, EndoCAM, CD31,
Immunogen	Mammalian cells - derived recombinant Human CD31 / PECAM1 (Met1-Lys601).
Accession	P16284
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.