

Product Information	
Catalog No.	MA10174
Product name	Anti-Human PIK3CG Monoclonal Antibody (1A128)
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A128
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	PtdIns-3-kinase subunit gamma, Phosphatidylinositol 4, 5-bisphosphate 3-kinase catalytic subunit gamma isoform, PtdIns-3-kinase subunit p110-gamma, PI3Kgamma, p120-PI3K, Phosphatidylinositol 4, 5-bisphosphate 3-kinase 110 kDa catalytic subunit gamma, PI3-kinase subunit gamma, PI3K-gamma, Phosphoinositide-3-kinase catalytic gamma polypeptide, Serine/threonine protein kinase PIK3CG, p110gamma, PIK3CG,
Immunogen	E. coli - derived recombinant Human PIK3CG (Trp624-Val718).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.