

| Product Information   |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | MA10180                                                                                                                                                      |
| Product name          | Anti-FITC/Fluorescein isothiocyanate Monoclonal Antibody (1A183)                                                                                             |
| Applications          | ELISA, IHC, WB                                                                                                                                               |
| Species reactivity    | General                                                                                                                                                      |
| Host species          | Mouse                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                   |
| Isotype               | IgG1                                                                                                                                                         |
| Clone ID              | 1A183                                                                                                                                                        |
| Form                  | Liquid                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                | FITC, Fluorescein isothiocyanate, 3326-32-7                                                                                                                  |
| Immunogen             | FITC conjugated with Carrier protein.                                                                                                                        |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Note                  | For research use only.                                                                                                                                       |