

Product Information	
Catalog No.	MA10181
Product name	Anti-Chicken OVA Monoclonal Antibody (1A178)
Applications	ELISA, IHC, WB
Species reactivity	Chicken
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A178
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	OVA, Ovalbumin, Allergen Gal d II, Egg albumin, Plakalbumin, Gal d 2
Immunogen	Natural OVA.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.