

Product Information	
Catalog No.	MA10196
Product name	Anti-Human GAPDH Monoclonal Antibody (1A200)
Applications	ELISA, WB
Species reactivity	Human, Mouse, Rat
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	1A200
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Protein A/G purified from cell culture supernatant.
Target	GAPD, Peptidyl-cysteine S-nitrosylase GAPDH, GAPDH, Glyceraldehyde-3-phosphate dehydrogenase
Immunogen	E. coli - derived recombinant Human GAPDH protein (Met1-Glu335).
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.