

Product Information	
Catalog No.	MA10227
Product name	Anti-Mouse TMEM87A Monoclonal Antibody (1A231)
Applications	ELISA, WB
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	1A231
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Transmembrane protein 87A, Elkin1, TMEM87A, PSEC0094
Immunogen	Recombinant Mouse TMEM87A protein.
Accession	Q8BXN9
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.