

Product Information	
Catalog No.	MA10229
Product name	Anti-SMCC Crosslinker Monoclonal Antibody (1A233)
Applications	ELISA
Species reactivity	General
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Clone ID	1A233
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	SMCC Crosslinker, Cys-conjugated, Succinimidyl-4-(N-maleimidomethyl)cyclohexane-1-carboxylate
Immunogen	Peptide and KLH complexes conjugated with SMCC Crosslinker.
Accession	CAS: 64987-85-5
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.