


Product Information	
Catalog No.	MA10242
Product name	Anti-HBV HBcAg Monoclonal Antibody (1A246)
Applications	ELISA
Species reactivity	Hepatitis B virus (HBV)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A246
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Capsid protein, Core antigen, Core protein, HBcAg, p21.5, C
Immunogen	HBV HBcAg Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human total PSA/tPSA Monoclonal Antibody (1A273)

SDS-PAGE