


Product Information	
Catalog No.	MA10249
Product name	Anti-Clostridium Difficile Toxin B/tcdB Monoclonal Antibody (1A253)
Applications	GICA, IFA
Species reactivity	Clostridioides difficile (Peptoclostridium difficile)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A253
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Toxin B, 3.4.22.-, Glucosyltransferase TcdB, 2.4.1.-, tcdB, toxB
Immunogen	Clostridium Difficile Toxin B/tcdB Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human MUC16/CA125 Monoclonal Antibody (1A280)

SDS-PAGE