


Product Information	
Catalog No.	MA10255
Product name	Anti-Influenza A virus/Flu-A Monoclonal Antibody (1A259)
Applications	GICA
Species reactivity	Influenza A virus (Flu-A)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A259
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Immunogen	Influenza A virus/Flu-A Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human HE4 Monoclonal Antibody (1A286)

SDS-PAGE