


Product Information	
Catalog No.	MA10257
Product name	Anti-Influenza B virus/Flu-B virus Monoclonal Antibody (1A261)
Applications	GICA
Species reactivity	Influenza B virus (Flu-B)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A261
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Immunogen	Influenza B virus/Flu-B virus Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human PG I/PGA Monoclonal Antibody (1A288)

SDS-PAGE