


Product Information

Catalog No.	MA10272
Product name	Anti-Human free PSA/fPSA Monoclonal Antibody (1A276)
Applications	ELISA, IFA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A276
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Prostate-specific antigen, PSA, 3.4.21.77, Gamma-seminoprotein, Semin, Kallikrein-3, P-30 antigen, Semenogelase, KLK3, APS
Immunogen	Human free PSA/fPSA Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human ENO2/NSE Monoclonal Antibody (1A303)

SDS-PAGE