


Product Information

Catalog No.	MA10278
Product name	Anti-Human Serum Ferritin/SF Monoclonal Antibody (1A282)
Applications	ELISA, GICA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A282
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Ferritin heavy chain, Ferritin H subunit, 1.16.3.1, Cell proliferation-inducing gene 15 protein, Ferritin heavy chain, N-terminally processed, FTH1, FTH, FTHL6, OK/SW-cl.84, PIG15
Immunogen	Human Serum Ferritin/SF Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human TSH Monoclonal Antibody (1A309)

SDS-PAGE