


Product Information

Catalog No.	MA10283
Product name	Anti-Human PG I/PGA Monoclonal Antibody (1A287)
Applications	ELISA, TIA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A287
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Pepsin A-3, 3.4.23.1, Pepsinogen-3, PGA3
Immunogen	Human PG I/PGA Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human LTF/Lactotransferrin Monoclonal Antibody (1A314)

SDS-PAGE