


### Product Information

|                       |                                                                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.           | MA10295                                                                                                                                                      |
| Product name          | Anti-Human CA19-9/CA199 Monoclonal Antibody (1A299)                                                                                                          |
| Applications          | IFA                                                                                                                                                          |
| Species reactivity    | Human                                                                                                                                                        |
| Host species          | Mouse                                                                                                                                                        |
| Clonality             | Monoclonal                                                                                                                                                   |
| Isotype               | IgG                                                                                                                                                          |
| Clone ID              | 1A299                                                                                                                                                        |
| Form                  | Liquid                                                                                                                                                       |
| Storage buffer        | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Purification          | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                | CA19-9, CA199, Carbohydrate antigen199, Carbohydrate antigen 19-9                                                                                            |
| Immunogen             | Human CA19-9/CA199 Protein.                                                                                                                                  |
| Stability and Storage | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Note                  | For research use only.                                                                                                                                       |

### Image Data



SDS-PAGE for Anti-Human HP/Haptoglobin Monoclonal Antibody (1A326)

SDS-PAGE