


Product Information

Catalog No.	MA10305
Product name	Anti-Human TSH Monoclonal Antibody (1A309)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A309
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Thyrotropin subunit beta, Thyroid-stimulating hormone subunit beta, TSH-B, TSH-beta, Thyrotropin beta chain, Thyrotropin alfa, TSHB
Immunogen	Human TSH Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human MYO Monoclonal Antibody (1A336)

SDS-PAGE