


Product Information	
Catalog No.	MA10318
Product name	Anti-Human SAA Monoclonal Antibody (1A322)
Applications	IFA, TIA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A322
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Serum amyloid A-1 protein, SAA, Amyloid protein A, Amyloid fibril protein AA, Serum amyloid protein A(2-104), Serum amyloid protein A(3-104), Serum amyloid protein A(2-103), Serum amyloid protein A(2-102), Serum amyloid protein A(4-101), SAA1
Immunogen	Human SAA Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human TNNI3/cTNI Monoclonal Antibody (1A349)

SDS-PAGE