


Product Information	
Catalog No.	MA10340
Product name	Anti-Human NT-ProBNP Monoclonal Antibody (1A344)
Applications	CLIA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A344
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Natriuretic peptides B, Brain natriuretic factor prohormone, preproBNP, proBNP, Gamma-brain natriuretic peptide, Iso-ANP, NT-proBNP, NT-pro-BNP, NT-proBNP(1-76), proBNP(3-108), NPPB
Immunogen	Human NT-ProBNP Protein.
Accession	P16860
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human HCG Monoclonal Antibody (1A371)

SDS-PAGE