


Product Information

Catalog No.	MA10346
Product name	Anti-Human D-Dimer Monoclonal Antibody (1A350)
Applications	CLIA, GICA, IFA, TIA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	1A350
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Fibrinogen alpha chain, Fibrinopeptide A, Fibrinogen alpha chain, FGA, Fibrinogen beta chain, Fibrinopeptide B, Fibrinogen beta chain, FGB, Fibrinogen gamma chain, FGG, PRO2061
Immunogen	Human D-Dimer Protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human AMH Monoclonal Antibody (1A377)

SDS-PAGE