

Product Information	
Catalog No.	MA10429
Product name	Anti-Romiplostim Monoclonal Antibody (1A432)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A432
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Thrombopoietin, C-mpl ligand, ML, Megakaryocyte colony-stimulating factor, Megakaryocyte growth and development factor, MGDF, Myeloproliferative leukemia virus oncogene ligand, THPO, MGDF
Immunogen	Recombinant Romiplostim protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Alternative Names	Romiplate, Romiplostim, Nplate
Note	For research use only.