

Product Information	
Catalog No.	MA10441
Product name	Anti-Human RAB39B Monoclonal Antibody (1A444)
Applications	ELISA, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	1A444
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Purification	Protein A/G purified from cell culture supernatant.
Target	Ras-related protein Rab-39B, RAB39B
Immunogen	Recombinant Human RAB39B protein.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.