


Product Information	
Catalog No.	PA10007
Product name	Anti-PSMD9 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Bos taurus
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	26S proteasome regulatory subunit p27, 26S proteasome non-ATPase regulatory subunit 9, PSMD9
Immunogen	E. coli - derived recombinant Human PSMD9 (Ser2-Arg223).
Accession	O00233
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.

Image Data



Western blot analysis was performed using anti-PSMD9 polyclonal antibody at 1ug/ml on various samples.
 Lane 1: NIH3T3 cell lysate
 Lane 2: HepG2 cell lysate
 Lane 3: K562 cell lysate
 Lane 4: Lncap cell lysate
 Lane 5: Hela cell lysate

Western blot