

Product Information	
Catalog No.	PA10022
Product name	Anti-CCL22/MDC Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	MDC(1-69), SCYA22, Macrophage-derived chemokine, CCL22, Stimulated T-cell chemotactic protein 1, CC chemokine STCP-1, Small-inducible cytokine A22, C-C motif chemokine 22, MDC
Immunogen	E. coli - derived recombinant Human CCL22/MDC (Gly25-Gln93).
Accession	O00626
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.