

Product Information	
Catalog No.	PA10026
Product name	Anti-SOCS2 Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Bos taurus, Human, Mouse, Rat, Pig
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CIS-2, STAT-induced STAT inhibitor 2, SSI-2, SOCS2, STAT2, CIS2, Cytokine-inducible SH2 protein 2, SSI2, Suppressor of cytokine signaling 2, SOCS-2
Immunogen	E. coli - derived recombinant Human SOCS2 (Met1-Val198).
Accession	O14508
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.