

Product Information	
Catalog No.	PA10028
Product name	Anti-CXCL11/I-TAC Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	IP-9, I-TAC, C-X-C motif chemokine 11, ITAC, Small-inducible cytokine B11, SCYB11, Beta-R1, Interferon gamma-inducible protein 9, H174, SCYB9B, CXCL11, Interferon-inducible T-cell alpha chemoattractant
Immunogen	E. coli - derived recombinant Human CXCL11/I-TAC (Phe22-Phe94).
Accession	O14625
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.