

Product Information	
Catalog No.	PA10029
Product name	Anti-TOR1A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Rat, Cricetus cricetus, Macaca fascicularis
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Torsin-1A, TOR1A, Torsin family 1 member A, DQ2, Dystonia 1 protein, TORA, TA, Torsin ATPase-1A, DYT1
Immunogen	E. coli - derived recombinant Human TOR1A (Val33-Ser250).
Accession	O14656
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.