

Product Information	
Catalog No.	PA10100
Product name	Anti-CS Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Mouse, Amblyrhynchus cristatus, Bos taurus, Danio rerio, Chicken, Human, Iguana iguana, Kajikia audax, Katsuwonus pelamis, Macaca fascicularis, Mesocricetus auratus, Rat, Pig, Tetraodon nigroviridis, Thunnus albacares, Thunnus obesus, Xenopus laevis, Xenopus tropicalis, Xiphias gladius
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	CS, Citrate synthase, mitochondrial, Citrate (Si)-synthase
Immunogen	E. coli - derived recombinant Human CS (Lys103-Ala394).
Accession	O75390
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.