

Product Information	
Catalog No.	PA10119
Product name	Anti-KRT33A Polyclonal Antibody
Applications	ELISA, IHC, WB
Species reactivity	Human, Mouse, Ovis aries, Pan troglodytes
Host species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Purification	Purified by antigen affinity column.
Target	Hair keratin, type I Ha3-I, KRTHA3A, K33A, HKA3A, Keratin, type I cuticular Ha3-I, Keratin-33A, KRT33A, HHA3-I
Immunogen	E. coli - derived recombinant Human KRT33A (Lys57-Glu363).
Accession	O76009
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.
Product Usage Information	ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000
Note	For research use only.